

501 Waterford, C Wing, C.D.Barfiwala Marg(Juhu Lane), Andheri(W), Mumbai 400 058. INDIA

Email tcmktg@sarex.com URLwww.sarex.com

## PRODUCT LITERATURE

# SUPERCONE-2100 (CONC)

## Product Information

Supercone-2100 (Conc) is a concentrated cold water soluble organic softener for cotton yarns, acrylic, hosiery, wool, silk and garments. Unlike conventional cationic softeners, Supercone-2100 (Conc) finished fabric shows minimum yellowing and shade change, hence it is ideal for finishing of whites and light dyed material.

## Key Features & Benefits

Key Features Benefits
----Concentrated organic Economical

softener

Minimum yellowing Suitable for finishing of whites

and pale shades

## General Characteristics

Physical appearance : Off white to creamish to light yellowish

flakes
: Cationic

Ionic nature : Cationic pH of 1% solution : 4 +/- 1

Miscibility : Dispersible with water

Compatibility : Compatible with cationic and nonionic

products

Stability : Stable to dilute acids and unstable to

dilute alkalies

## Application

Padding process Exhaust process

Dosage\* : 4-7 g/l Dosage\* : 0.4-1%
Pick-up : 65-70% Bath pH : 5.0-6.0
Bath pH : 5.0-6.0 Bath Temp. : 30-40 deg.C
Drying : 130-160 deg.C Time : 20-30 mins.

\* Dosages recommended are based on the concentrated product's strength.

501 Waterford, C Wing, C.D.Barfiwala Marg(Juhu Lane), Andheri(W), Mumbai 400 058. INDIA

Email tcmktg@sarex.com URLwww.sarex.com

## PRODUCT LITERATURE

# SUPERCONE-2100 (CONC)

Dosages of diluted product to be calculated on the basis of solid content.

### Instructions For Dilution

Procedure to prepare 10% diluted product of 100% Supercone-2100 (Conc)

#### Method-I

#### \_\_\_\_\_

- Take 89 parts of DM water. Adjust the pH of water to 4.5-5.5
- To this, add 11 parts of Supercone-2100 (Conc). Soak for 8-12 hrs.
- After 8-12 hrs, stir constantly for 20-30 mins at 300-400 rpm to obtain homogeneous product

#### Method-II

#### \_\_\_\_\_

- Take 89 parts of DM water. Adjust the pH of water to 4.5-5.5
- To this, gradually add 11 parts of Supercone-2100 (Conc) with constant stirring at 300-400 rpm
- Stir for 30-45 mins to obtain homogeneous product and allow standing for overnight to completely dissolve all the particles.

### Preservation:

#### -----

It is recommended to add 0.1% Bronopol (2-Bromo,2-nitropropane, 1,3-diol) or Methyl Paraben (methyl para hydroxy benzoate) as an antimicrobial agent to avoid fungal/microbial growth.

Note: Before use, stir for 5 min.

It is advisable to utilize the diluted product within 2 months. Stir well before use.

#### Precautions

Storage

: Store in cool, ventilated shed away from heat and direct sunlight. Storage temperature should not exceed 35 deg C. Close lids firmly to avoid contact with air and moisture.



501 Waterford, C Wing, C.D.Barfiwala Marg(Juhu Lane), Andheri(W),Mumbai 400 058. INDIA

Email tcmktg@sarex.com URLwww.sarex.com

# **PRODUCT LITERATURE**

# SUPERCONE-2100 (CONC)

DOLDICOND		(00110)		
01-1 C T ' C-	0	C		
Shelf Life	: 9 months stored	from the date of under controlled	manufacturing, conditions.	lİ
	202200	diido1 00110101100	•	

This is a computer generated report. Hence not signed.